

## FlexWell™ Silicone Isolator, Seal Strips, and Tray

Catalog nos. F24781, F24782, F24783


**Table 1.** FlexWell™ products.

Product	Dimensions	Volume	Depth	Quantity per Package
FlexWell™ silicone isolator with adhesive, 16-wells				
F24781	25 mm × 75 mm 6.8 mm × 6.8 mm (well dimensions)	15–150 µL	3.3 mm	10
FlexWell™ tray and cover				
F24782	85.4 mm × 127.5 mm	NA	NA	1
FlexWell™ seal strips and applicator				
F24783	25 mm × 86 mm	NA	NA	50 strips and 1 applicator
NA: not applicable.				

### Introduction

FlexWell™ silicone isolator is a disposable gasket that fits onto a standard microscope slide to form wells to isolate up to 16 specimens per slide. The slide with FlexWell™ silicone isolator fits into a reusable FlexWell™ tray and cover. Four FlexWell™ silicone isolators fit into a tray, producing a modular plate with a standard microtiter plate footprint and well spacing allowing automated robotic processing. The 9 mm center to center well spacing facilitates use with multichannel pipettors.

The FlexWell™ seal strips and applicator are removable adhesive covers for FlexWell™ silicone isolators and is used to seal for mixing and to prevent evaporation.



**Figure 1.** FlexWell™ silicone isolators fitted in a FlexWell™ tray.

## Product List Current prices may be obtained from our website or from our Customer Service Department.

Cat. no.	Product Name	Unit Size
F24781	FlexWell™ silicone isolator with adhesive, disposable 16 well gasket *set of 10*	1 set
F24782	FlexWell™ tray and cover, reusable, holds up to 4 slides *includes one tray and cover*	1 set
F24783	FlexWell™ seal strips and applicator, removable adhesive covers *set of 50 seal-strips and one applicator*	1 set

## Contact Information

### Molecular Probes, Inc.

29851 Willow Creek Road  
Eugene, OR 97402  
Phone: (541) 465-8300  
Fax: (541) 335-0504

### Customer Service:

6:00 am to 4:30 pm (Pacific Time)  
Phone: (541) 335-0338  
Fax: (541) 335-0305  
probesorder@invitrogen.com

### Toll-Free Ordering for USA:

Order Phone: (800) 438-2209  
Order Fax: (800) 438-0228

### Technical Service:

8:00 am to 4:00 pm (Pacific Time)  
Phone: (541) 335-0353  
Toll-Free (800) 438-2209  
Fax: (541) 335-0238  
probestech@invitrogen.com

### Invitrogen European Headquarters

Invitrogen, Ltd.  
3 Fountain Drive  
Inchinnan Business Park  
Paisley PA4 9RF, UK  
Phone: +44 (0) 141 814 6100  
Fax: +44 (0) 141 814 6260  
Email: euroinfo@invitrogen.com  
Technical Services: eurotech@invitrogen.com

**For country-specific contact information, visit [www.invitrogen.com](http://www.invitrogen.com).**

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

### Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2009, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.