

Reactive Oxygen Species (ROS) Detection Reagents

Table 1. Contents and storage information.

Material	Amount	Storage	Stability
Reactive Oxygen Species (ROS) Detection reagents	See label	<ul style="list-style-type: none"> • ≤ -20°C • Desiccate • Protect from light • Protect from air * • Avoid freeze-thaw cycles 	When stored as directed, this product is stable for at least 3 months.

* These products are air sensitive and should be stored under dry argon or nitrogen.

Approximate Fluorescence Excitation and Emission, in nm : 492–495/517–527, in nm.

Introduction

Molecular Probes offers derivatives of reduced fluorescein and calcein as cell-permeant indicators for reactive oxygen species. Chemically reduced and acetylated forms of 2',7'-dichlorofluorescein (DCF) and calcein are nonfluorescent until the acetate groups are removed by intracellular esterases and oxidation occurs within the cell^{1–4} (Figures 1 and 2). Esterase cleavage of the lipophilic blocking groups yields a charged form of the dye that is much better retained by cells than is the parent compound. Oxidation of these probes can be detected by monitoring the increase in fluorescence with a flow cytometer, fluorometer, microplate reader, or fluorescence microscope, using excitation sources and filters appropriate for fluorescein (FITC).^{5–8} Because the dyes are susceptible to photo-oxidation, low light conditions should be used for fluorescence microscopy applications whenever possible.⁹

The carboxy derivative of fluorescein, carboxy-H₂DCFDA (C400), carries additional negative charges that improve its retention compared to noncarboxylated forms. The fluorinated derivative, H₂DFFDA (C13293), exhibits improved photostability compared to chlorinated fluorescein derivatives. Derivatives with a thiol-reactive chloromethyl group (C6827) or an amine-reactive succinimidyl ester group (D2935) allow for covalent binding to intracellular components, permitting even longer retention within the cell. H₂DCFDA, SE (D2935) may also provide extracellular labeling.

Dihydrocalcein, AM is freely permeant to cell membranes and is oxidized to green-fluorescent calcein, which has excellent retention properties. However, unlike the fluorescein derivatives, calcein fluorescence can be quenched by Fe³⁺, Co²⁺, Ni²⁺, and Cu²⁺ at pH 7. Other reduced fluorescein ROS indicators available from Molecular Probes include aminophenyl fluorescein (A36003) and hydroxyphenyl fluorescein (H36004).




Figure 1. Absorption spectrum of reduced dye (---) and absorption/emission spectra of oxidized dye (—).

Guidelines for Use

Materials Required but Not Provided

- DMSO
- DMF
- Ethanol
- Buffer, such as PBS, HBSS, or HEPES

Preparing the Stock Solutions

Do not dissolve products until immediately prior to use. Acetylated dyes should be reconstituted only as required using high quality anhydrous dimethylsulfoxide (DMSO), dimethylformamide (DMF), or 100% ethanol. Working solutions should be freshly prepared; discard excess diluted probe at the end of the work session. These probes oxidize more readily in solution, and the presence of moisture will facilitate the decomposition of the dye.

General Considerations

Protocols for loading a particular dye and/or cell type can be found in the literature. The AM or acetate ester concentration should be kept as low as possible to reduce potential artifacts from overloading, including incomplete hydrolysis, compartmentalization, and the toxic effects of hydrolysis by-products (e.g., formaldehyde, acetic acid). Loading may be done at the optimal temperature for cell growth; however, some investigators have reported a greater degree of compartmentalization at physiological temperatures than is observed at room temperature.




Figure 2. The structures of (a) H₂DCFDA, (b) carboxy-H₂DCFDA, (c) chloromethyl-H₂DCFDA, (d) deacetylated H₂DCF, and (e) the deacetylated, oxidized product, DCF.

To determine the optimal dye-loading concentration, temperature, incubation, recovery, and retention time for your specific application, we recommend using a cell-permeant oxidized form of the ROS indicator such as carboxy-DCFDA (C369), fluorescein diacetate (F1303), or calcein, AM (C1430).

Extracellular hydrolysis of the AM and acetate esters can be minimized by using a loading buffer that is free of primary and secondary amines (such as PBS). Cells should be maintained in medium that is free of phenol red and other colorimetric dyes prior to and throughout the assay and washed with dye-free buffer or media after loading.

Note: The time the dye is retained within the cells must cover the duration of the experiment. Cells vary in their response to dyes; some dyes are lost by passive diffusion over time, or more quickly by active efflux. If the assay is to be carried over an extended period of time, dye-retention times need to be determined beforehand. This should also be done for derivatives that allow for covalent attachment, since cells vary in the rates at which they turn over various cellular components, and in their viability subsequent to covalent modification.

Protocols The following protocols provide general guidelines derived from various publications, and should be modified for the particular application and sensitivity required.

- 1.1 Shortly before performing the experiments, reconstitute the ROS indicator to make a concentrated stock solution. Keep tightly sealed until ready to use.
- 1.2 Remove cells from growth media via centrifugation or pipetting. Resuspend cells in pre-warmed buffer (PBS, HBSS, HEPES, or other simple physiological buffer) containing the probe to provide a final working concentration of ~1–10 μM dye. The optimal working concentration for your application must be empirically determined.
- 1.3 Incubate at the optimal temperature for the cells. Generally, a loading time of 5–60 minutes is sufficient.
- 1.4 Remove the loading buffer; return the cells to prewarmed growth medium and incubate at the optimal temperature. For derivatives with acetoxymethyl ester (AM) and/or diacetate groups, allow a short recovery time for cellular esterases to hydrolyze the AM or acetate groups and render the dye responsive to oxidation. The optimal recovery time can vary widely, as some cell types normally exhibit low levels of esterase activity.¹⁰
- 1.5 Determine the baseline fluorescence intensity of a sample of the loaded cells prior to exposing the cells to experimental inducements. If using the succinimidyl ester derivative (C2935), any extracellularly bound dye can be quenched using Trypan Blue (~0.0025% w/v) in order to better distinguish the signal from the intracellular ROS response.
- 1.6 Negative controls should be assessed as follows:
 - a) Examine unstained cells for autofluorescence in the green emission range.
 - b) For flow cytometry, ascertain that the forward and side scatter of cells is unchanged after dye-loading and treatment. Changes in cell dimensions may be related to blebbing or shrinkage resulting from handling or a toxic response.
 - c) Examine the fluorescence of cell-free mixtures of dye and buffer/media with and without the inducer/stimulant. In the absence of extracellular esterases and other oxidative enzymes, the ROS indicator should exhibit a gradual increase in fluorescence over time, which may be related to spontaneous hydrolysis, atmospheric oxidation, and/or light-induced oxidation.^{11–13}
 - d) Examine the fluorescence of untreated loaded cells that have been maintained in growth medium or simple buffer. In healthy cells, oxygen radicals are eliminated by cellular enzymes and/or natural antioxidants. Following the dye-loading recovery period, healthy cells should exhibit a low level of fluorescence that is relatively stable for the duration of the experiment;

however, a gradual increase (due to auto-oxidation) or decrease (due to loss of dye from cells or photobleaching) in fluorescence may be observed. In the absence of any stimulus or inducement, a burst of fluorescence in healthy, untreated cells could indicate progress to cell death or some other oxidative event.

1.7 To create positive controls, oxidative activity may be stimulated with:

a) the tumor promoter 4 β -phorbol 12-myristate 13-acetate (PMA; stock solution 1 mM in DMF; working concentration 100 pM to 10 μ M)

b) the bacterial chemotactic peptide N-formyl-L-methionyl-L-leucyl-L-phenylalanine (fMLP; stock solution 1 mM in DMF; working concentration 1–10 μ M)

c) H₂O₂ or *tert*-butyl hydroperoxide (TBHP) to a final concentration of ~100 μ M (increase or decrease based on the sensitivity and response of the cells).

1.8 To assure that an added drug or compound will not quench the dye, examine the absorbance spectrum of the compound and determine that the absorbance peak does not overlap with either the excitation or emission peak of the oxidized dye. Alternatively, you can mix a solution of the drug/compound with carboxy-DCF (C368), fluorescein (F1300), or calcein (C481), then compare the fluorescence intensity to a control solution of the dye, or subject a culture loaded with the cell-permeant, oxidized form of the dye to the drug/compound and compare to dye-loaded cells untreated.

References

1. Anal Biochem 11, 6 (1965); **2.** Anal Biochem 11, 1 (1965); **3.** Anal Biochem 134, 111 (1983); **4.** Cell Biol Int 24, 757 (2000); **5.** J Immunol Methods 156, 39 (1992); **6.** J Immunol Methods 178, 89 (1995); **7.** Toxicol in Vitro 11, 531 (1997); **8.** J Histochem Cytochem 50, 289 (2002); **9.** Biochem Biophys Res Commun 304, 619 (2003); **10.** J Leukoc Biol 43, 304 (1988); **11.** FEBS Lett 511, 21 (2002); **12.** Free Radic Biol Med 33, 938 (2002); **13.** Free Radic Biol Med 30, 1418 (2001).

Product List

Current prices may be obtained from our website or from our Customer Service Department.

Cat #	Product Name	Unit Size
C13293	5-(and-6)-carboxy-2',7'-difluorodihydrofluorescein diacetate (carboxy-H ₂ DFFDA) *mixed isomers*	5 mg
C400	5-(and-6)-carboxy-2',7'-dichlorodihydrofluorescein diacetate (carboxy-H ₂ DCFDA) *mixed isomers*	25 mg
C6827	5-(and-6)-chloromethyl-2',7'-dichlorodihydrofluorescein diacetate, acetyl ester (CM-H ₂ DCFDA) *mixed isomers* *special packaging*	20 x 50 μ g
D23805	dihydrocalcein, AM *special packaging*	20 x 50 μ g
D2935	2',7'-dichlorodihydrofluorescein diacetate, succinimidyl ester (OxyBURST® Green H ₂ DCFDA, SE)	5 mg
D399	2',7'-dichlorodihydrofluorescein diacetate (2',7'-dichlorofluorescein diacetate; H ₂ DCFDA)	100 mg
C2938	6-carboxy-2',7'-dichlorodihydrofluorescein diacetate, di(acetoxymethyl ester)	5 mg

Contact Information

Molecular Probes, Inc.
29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:
6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:
Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:
8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probetech@invitrogen.com

Invitrogen European Headquarters
Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Assistance Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2006, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.