

Vps10p (yeast) monoclonal antibody

CATALOG #:	A21274
COMPONENTS:	250 µg monoclonal antibody
APPLICATIONS:	Western blotting
CLONE ID OF MONOCLONAL ANTIBODY (mAb):	18C8BC7AD10
SPECIES CROSS-REACTIVITY:	<i>S. cerevisiae</i>
HOST SPECIES AND ISOTYPE:	Mouse IgG2a, k
IMMUNOGEN:	Yeast Vps10p
CONCENTRATION:	0.5 mg/mL in HEPES-Buffered Saline (HBS) with 0.02% azide as a preservative
SUGGESTED WORKING CONCENTRATION:	4.0 µg/mL for Western blotting
mAb PURITY:	Near homogeneity as judged by SDS-PAGE. The antibody was produced <i>in vitro</i> using hybridomas grown in serum-free medium, and then purified by biochemical fractionation.
STORAGE CONDITIONS:	Store at 4°C. Do not freeze.

For research use only. Not intended for human or animal therapeutic or diagnostic use.

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

www.invitrogen.com