

# Complex IV subunit VIb monoclonal antibody

<b>CATALOG #:</b>	A21366
<b>COMPONENTS:</b>	100 µg monoclonal antibody
<b>APPLICATIONS:</b>	Western blotting
<b>CLONE ID OF MONOCLONAL ANTIBODY (mAb):</b>	3F9D3D11AF6
<b>SPECIES CROSS-REACTIVITY:</b>	human, bovine, and rat
<b>HOST SPECIES AND ISOTYPE:</b>	Mouse IgG1, k
<b>IMMUNOGEN:</b>	Bovine Complex IV (native holoenzyme protein, purified from heart)
<b>CONCENTRATION:</b>	1 mg/mL in Hepes-Buffered Saline (HBS) with 0.02% azide as a preservative
<b>SUGGESTED WORKING CONCENTRATION:</b>	1 µg/mL for Western blotting
<b>mAb PURITY:</b>	Near homogeneity as judged by SDS-PAGE. The antibody was produced <i>in vitro</i> using hybridomas grown in serum-free medium, and then purified by biochemical fractionation.
<b>STORAGE CONDITIONS:</b>	Store at 4°C. Do not freeze.

## BACKGROUND:

Complex IV, also called cytochrome *c* oxidoreductase or cytochrome *c* oxidase, is a complex of 13 different subunits, three of which (I, II and III) are encoded on mitochondrial DNA and the remainder in the nuclear DNA. The complex contains two heme groups (*a* and *a*<sub>3</sub>) and two copper atoms as prosthetic groups. Genetic alterations of this enzyme complex are a common cause of OXPHOS diseases and the enzyme is altered in patients with Alzheimer's disease. Also, there are reports of reduced amounts of this complex in hypoxic cancer cells.

**For research use only. Not intended for human or animal therapeutic or diagnostic use.**

**Important Licensing Information** - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, [www.invitrogen.com](http://www.invitrogen.com)). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

[www.invitrogen.com](http://www.invitrogen.com)

Life Technologies Corporation • 5791 Van Allen Way • Carlsbad • CA 92008 • Tel: 800.955.6288 • E-mail: [techsupport@invitrogen.com](mailto:techsupport@invitrogen.com)

MP21366

MAN0003198

Rev. date: 12 July 2010