

69 kDa V-ATPase monoclonal antibody

CATALOG #:	A6422
COMPONENTS:	250 µg monoclonal antibody
APPLICATIONS:	Western blotting and Immunocytochemistry (heat-induced antigen-retrieval improves signal)
CLONE ID OF MONOCLONAL ANTIBODY (mAb):	8B1F3
SPECIES CROSS-REACTIVITY:	<i>S. cerevisiae</i>
HOST SPECIES AND ISOTYPE:	Mouse IgG2a, k
IMMUNOGEN:	Purified <i>S. cerevisiae</i> vacuole membranes
CONCENTRATION:	1 mg/mL in HEPES-Buffered Saline (HBS) with 0.02% azide as a preservative
SUGGESTED WORKING CONCENTRATION:	0.1 µg/mL for Western blotting. 20 µg/mL for Immunocytochemistry
mAb PURITY:	Near homogeneity as judged by SDS-PAGE. The antibody was produced <i>in vitro</i> using hybridomas grown in serum-free medium, and then purified by biochemical fractionation.
STORAGE CONDITIONS:	Store at 4°C. Do not freeze.
COUNTRY OF ORIGIN:	USA

BACKGROUND:

Antibody binding specificity was determined by Particle Concentration Fluorescence Immunoassay (PCFIA) using the 69 kDa vacuolar H⁺-ATPase (V-ATPase) subunit, by Western blot immunoassay using protein extracts from yeast, and by indirect immunofluorescence of fixed yeast cells.

The antibody can be used to detect the yeast 69 kDa V-ATPase subunit (the product of VMA1/TFP1 gene) on Western blots or to immunolocalize the V-ATPase enzyme complex in fixed yeast cells. The antibody can also be used to precipitate the native yeast 6/6 V-ATPase complex.

For research use only. Not intended for human or animal therapeutic or diagnostic use.

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

www.invitrogen.com

Life Technologies Corporation • 5791 Van Allen Way • Carlsbad • CA 92008 • Tel: 800.955.6288 • E-mail: techsupport@invitrogen.com

MP64220

MAN0002728

Rev. date: 01 June 2010