

Change in Packaging for Products:

P-3584 POPO™-3 iodide (534/570)

N-7565 Nucleic Acid Stains Dimer Sampler Kit (N-7565)

Dear Customer:

Molecular Probes has recently changed the solvent that we use for product POPO™-3 iodide (534/570) (P-3584). The product is now offered as a 1 mM solution in *N,N*-dimethylformamide (DMF); formerly it was provided as a 1 mM solution in dimethylsulfoxide (DMSO). This packaging change should not affect the performance of the dye and should improve its stability upon storage. This change also applies to the sample size of POPO-3 dye, Component D, in our Nucleic Acid Stains Dimer Sampler Kit (N-7565).

We apologize for any inconvenience that this change may have caused. If you have any questions about this product, please contact our Technical Assistance Department.